



# DELIVER CARE SAFELY.

**Medication Delivery and Disposal**  
Ensure end-to-end safe delivery and disposal of medications

# Needlestick injuries don't need to happen.

Exposure to bloodborne pathogens. Loss of employee time. Cost of treatment for infected staff. **Preventing and treating sharps-related injuries is more challenging than ever.**

What do you need to succeed? A hands-on medical supplies partner who shares your passion for providing safe and efficient care.

Our comprehensive product line can help you provide end-to-end safe delivery and disposal of medications. We supply a full range of products from needles and syringes to IV catheters and sharps containers.

## The numbers tell the story

**385K**

Sharps-related injuries that occur annually among healthcare workers in hospitals.<sup>1</sup>

**31%**

Percent of sharps injuries to healthcare workers caused by disposable syringes.<sup>1</sup>

**\$3,042**

Estimated cost per needlestick (includes: lab fees for testing exposed employees, associated labor and counseling, and post-exposure follow-ups).<sup>2</sup>

## SureSite™ IV Safety Catheter



DYNSCS22100

### SureSite™ IV Slide Safety Catheters

- » Self-sheathing needle shield with a clicking lock safety mechanism
- » Flashback chamber provides clear visual confirmation
- » Needle bevel designed for smooth insertion and greater patient comfort

Item No.	Description	Pkg
DYNSCS18114	SureSite IV Slide Catheter, 18 G x 1¼"	200/cs
DYNSCS20114	SureSite IV Slide Catheter, 20 G x 1¼"	200/cs
DYNSCS20100	SureSite IV Slide Catheter, 20 G x 1"	200/cs
DYNSCS22100	SureSite IV Slide Catheter, 22 G x 1"	200/cs
DYNSCS24340	SureSite IV Slide Catheter, 24 G x ¾"	200/cs



DYNCA22100

### Suresite™ IV Auto Safety Catheters

- » Fully encapsulated needle shield and push button safety activation
- » Flashback chamber provides clear visual confirmation
- » Needle bevel designed for smooth insertion and greater patient comfort

Item No.	Description	Pkg
DYNCA18116	SureSite IV Auto Catheter, 18 G x 1.16"	200/cs
DYNCA18188	SureSite IV Auto Catheter, 18 G x 1.88"	200/cs
DYNCA20100	SureSite IV Auto Catheter, 20 G x 1"	200/cs
DYNCA20116	SureSite IV Auto Catheter, 20 G x 1.16"	200/cs
DYNCA20188	SureSite IV Auto Catheter, 20 G x 1.88"	200/cs
DYNCA22100	SureSite IV Auto Catheter, 22 G x 1"	200/cs
DYNCA24075	SureSite IV Auto Catheter, 24 G x 0.75"	200/cs

## IV Administration



DYNDTC5081

### IV Extension Sets

- » Luer lock connectors allow fluid to flow from the bag to the patient; cap pops off easily with the push of a thumb
- » Slide clamps eliminate air bubbles and temporarily stops fluid flow to avoid resetting the roller clamp flow rate
- » Non-bonded needle-free valve

Item No.	Description	Pkg
DYNDTC5081	8" Set with Nonbonded Needle-Free Valve, 0.3 mL	100/cs
DYNDTC5077	6.5" Set with Nonbonded Needle-Free Valve, 0.7 mL	100/cs
DYNDTC5082	8" Set with Needle-Free T Connector, Small Bore, 0.4 mL	100/cs
DYNDTN1002	6" Set, 0.2 mL, M/F Luer Lock Connectors, Side Clamps	100/cs
DYNDTN0007	60" Set, Mini Bore, 0.6 mL	100/cs
DYNDTN1603	60" Set, Small Bore, 2.0 mL	100/cs



DYNDTN1545

### IV Administration Sets

- » Universal drip chamber features wide finger flanges and a very sharp spike for quickly puncturing the IV bag
- » Non-removable slide clamp eliminates air bubbles and temporarily stop fluid flow to avoid resetting the roller clamp flow rate
- » Rotating male luer lock

Item No.	Description	Pkg
DYNDTN1529	15DRP, 1 Needle-Free Port, 92", PV 21 mL	50/cs
DYNDTN1530	15DRP, 1 Injection Port, 92", PV 21 mL	50/cs
DYNDTN1512	15DRP, 2 Needle-Free Ports, 96", PV 21 mL	50/cs
DYNDTN0512	60DRP, 2 Needle-Free Ports, 96", PV 22 mL	50/cs
DYNDTN1545	15DRP, 3 Needle-Free Ports, 109", PV 23 mL	50/cs
DYNDTN1540A	Secondary, 40" with Hanger, PV 13 mL	50/cs
DYNDTC2000A	SmartSite Needle-Free Valve	100/cs

## Needles and Syringes



SYRS103255

### Syringes with Safety Needles

- » Preassembled needle-and-syringe combinations with features to help ensure safety for nurses and patients

Item No.	Description	Pkg
SYRS103227	Syringe with Safety Needle, 22 G x 1.5", 3 mL	1200/cs
SYRS103255	Syringe with Safety Needle, 25 G x 1", 3 mL	1200/cs
SYRS103253	Syringe with Safety Needle, 25 G x 5/8", 3 mL	1200/cs
SYRS101272F	Syringe with Safety Needle, 27 G x 0.5", 1 mL	1200/cs



SYRS100253

### Safety Hypodermic Needles

- » Advanced and reliable features ensure safety for nurses and patients

Item No.	Description	Pkg
SYRS100185	Safety Hypodermic Needle, 18 G x 1"	800/cs
SYRS100235	Safety Hypodermic Needle, 23 G x 1"	800/cs
SYRS100255	Safety Hypodermic Needle, 25 G x 1"	800/cs
SYRS100253	Safety Hypodermic Needle, 25 G x 5/8"	800/cs



SYRSI101292

### Insulin Safety Syringes

- » Protective shield temporarily covers needle; after injection, the shield is permanently locked in place by pulling forward and twisting
- » Compliant with OSHA Bloodborne Pathogens Standard
- » Low dead-space design

Item No.	Description	Pkg
SYRSI101292	Insulin Safety Syringe, 29 G x 0.5", 1 mL	500/cs
SYRSI193292	Insulin Safety Syringe, 29 G x 0.5", 0.5 mL	500/cs



SYR110022

### Standard and Blunt Fill Needles and Combinations

- » Blunt fill needles with and without filter designed to draw from vials and prevent coring
- » High quality standard needles and combos for general hypodermic purposes

Item No.	Description	Pkg
SYR110022	Blunt Needle, 18 G x 1.5"	1000/cs
SYR110022F	Blunt Needle with Filter, 18 G x 1.5"	1000/cs
SYR100187	Standard Hypodermic Needle, 18 G x 1.5"	1000/cs
SYR103205	Standard Hypodermic Syringe with Needle, LL, 20 G x 1", 3 mL	800/cs

### Allergy Syringe Trays

- » Smooth-gliding rubber stopper promotes accurate medication measurement
- » Low dead-space minimizes allergen waste
- » Easy-to-read number scale helps you prepare precise dose
- » Sterile packaged 25 per tray for easy accessibility; trays stack for ease and convenience

SYRA101272



Item No.	Description	Pkg
SYRA101272	Allergy Tray, 27 G x 0.5", 1 mL	40 trays/cs
SYRA101261	Allergy Tray, 26 G x 3/8", 1 mL	40 trays/cs



SYR110010

### Luer Lock Syringes

- » Manufactured to precise tolerance
- » Available in 1-60 mL
- » Sterile and individually packaged

Item No.	Description	Pkg
SYR103010	Luer Lock Sterile Syringe, 3 mL	800/cs
SYR105010	Luer Lock Sterile Syringe, 5 mL	400/cs
SYR110010	Luer Lock Sterile Syringe, 10 mL	400/cs
SYR120010	Luer Lock Sterile Syringe, 20 mL	160/cs
SYR160010	Luer Lock Sterile Syringe, 60 mL	160/cs

## Blood Collection

### OnGuard™ Blood Collection Sets

- » Large anti-skid wings for enhanced control
- » Sturdy non-slip grip to engage the safety feature
- » Compliant with OSHA Bloodborne Pathogens Standard

Item No.	Description	Pkg
DYNDTNBCS21G	Blood Collection Set, 21 G x 0.75"	200/cs
DYNDTNBCS23G	Blood Collection Set, 23 G x 0.75"	200/cs
DYNDTNBCS25G	Blood Collection Set, 25 G x 0.75"	200/cs



DYNDTNBCS23G

### OnGuard™ Blood Collection Holders

- » Helps your institution comply with OSHA regulations
- » Grooves for easy handling
- » Compatible with most blood collection sets, needles and tubes

Item No.	Description	Pkg
DYNDBCHLDR	Blood Collection Tube Holder	1000/cs



DYNDBCHLDR

## Sharps Containers

### Biohazard Patient Room Sharps Containers

- » Counter-balanced lid creates a torturous path to facilitate safe one-handed disposal
- » Horizontal drop design for maximum capacity to fill line
- » Each container locks for final disposal
- » Available in wall mount and freestanding designs with multiple bracket and cabinet options
- » Easy-to-see fill levels

Item No.	Description	Pkg
MDS705153	Sharps Container, 5 Qt, Translucent Red, Counter-Balance Lid	20/cs
MDS705253	Sharps Container, 5.4 Qt, Translucent Red, Counter-Balance Lid	30/cs
MDS707953	Locking Cabinet for MDS705153	5/cs
MDS705203	Sharps Container, 3 Gal, Red	12/cs



MDS705153 (Shown with cabinet MDS707953)



## Phlebotomy & Lab Containers

- » Designed to fit into blood drawing trays for point-of-use disposal
- » Square base for added stability
- » Available with screw top or dual opening lids
- » Brackets available to mount containers to wall and free up counter space are available separately

Item No.	Description	Pkg
MDS705110	1 Qt, Red, Phlebotomy Sharps Container	100/cs
MDS705115	1.5 Qt, Red, Transportable Container	20/cs
MDS705202	2 Gal, Red, Sliding Lid	20/cs



MDS705202

MDS705110

## Large Sharps Containers

- » These containers are designed for use in restricted access areas
- » Variety of lid designs offer new to market options for all sizes of sharps
- » Integrated side handles make the container easy to carry
- » Rolling carts available for easy transport of large containers

Item No.	Description	Pkg
MDS705208HN	8 Gal, Red, Hinge Lid	10/cs
MDS705218HN	18 Gal, Red, Hinge Lid	5/cs
MDS705218SL	18 Gal, Red, Slide Lid	5/cs
MDS707CART	Rolling Trolley Cart for Large Biohazard Containers	1/cs



## Large PG-II Containers

- » Containers tested to comply with Department of Transportation (DOT) requirements to eliminate the need for secondary packaging upon disposals

Item No.	Description	Pkg
MDS706308	Large PG-II Containers, 8 Gal, Yellow	9/cs
MDS706318	18 Gal, Yellow, PGII Chemo Container, Port Lid	7/cs
MDS707318	18 Gal, White/Blue, PGII Pharmaceutical Waste, Port Lid	7/cs



MDS707318

MDS706318

**Learn more about our full line of products for medication delivery and disposal.**  
Contact your Medline Representative or go to [medline.com](http://medline.com).



**Medline Industries, Inc.** Three Lakes Drive, Northfield, IL 60093

**Medline United States**  
1-800-MEDLINE (633-5463)  
[medline.com](http://medline.com) | [info@medline.com](mailto:info@medline.com)

**Medline Canada**  
1-800-396-6996  
[medline.ca](http://medline.ca) | [canada@medline.com](mailto:canada@medline.com)

**Medline México**  
01-800-831-0898  
[medlinemexico.com](http://medlinemexico.com) | [mexico@medline.com](mailto:mexico@medline.com)

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

1. The National Institute for Occupational Safety and Health (NIOSH) website, <https://www.cdc.gov/niosh/stopsticks/sharpsinjuries.html>. Accessed November 11, 2018.  
2. Needlestick Injuries Cost U.S. Hospitals \$1B, Becker's, August 27, 2013. <https://www.beckershospitalreview.com/quality/needlestick-injuries-cost-u-s-hospitals-1b.html>. Accessed November 1, 2018. Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.  
©2018 Medline Industries, Inc. Medline is a registered trademark and OnGuard is a trademark of Medline Industries, Inc. MKT19W44328 / LIT530R / 2.5M / SG / 49